

## AP006: GETTING THE BEST OUT OF YOUR CAPACITIVE FORCE & DISPLACEMENT SENSORS

Razorbill instruments stress and strain cells use capacitors to provide a feedback signal of the force or displacement your sample is subject to. This app note covers the expected performance of the capacitors, and how to design your experiments get the best possible performance from them.

App note *AP003: Measuring Capacitance* covers suitable measuring devices and advice on connecting them to your cell, while app note *AP004: Cables and Heat load* contains advice on adding new cables to a cryostat where this is required. App note *AP007 Force or Displacement* covers the merits of the different types of sensors and how the value measured relates to the stress and strain in the sample.

# CONTENTS

- Introduction..... 3
  - Converting capacitance to force or displacement ..... 3
  - Methods for on-site (re)calibration ..... 4
- Factors affecting the calibration..... 5
  - Temperature ..... 5
  - Rate of change of temperature ..... 7
  - Exposure to atmosphere ..... 8
  - Hysteresis ..... 10
  - Effect of displacement on a force sensor ..... 11
  - Effect of force on a displacement sensor ..... 11
  - Other sources of capacitance change ..... 12
- Advice on experiments ..... 12
  - Finding the location of zero ..... 12
  - Converting capacitance to displacement or force, correcting for temperature and other factors. .... 13
- Troubleshooting ..... 14
  - Small steps in capacitance (less than 5 %)..... 14
  - Measured capacitance falls to zero ..... 15
  - Capacitor doesn't respond to changes in piezo voltage ..... 15
  - Hysteresis ..... 16
  - Locating a fault ..... 16
  - Helping us to help you ..... 16



# INTRODUCTION

Razorbill's feedback sensors rely on the variation in the capacitance of a parallel-plate capacitor as the plates are moved towards and away from each other. The design of the capacitor varies slightly, but they typically have an area of a few square mm and an initial spacing of some tens of microns, leading to a nominal capacitance of around a picofarad. To confirm the characteristics of the capacitor in your device, please refer to the datasheet and calibration curve for your specific cell.

As our smallest cells have working ranges of only a few microns, the change in capacitance during an experiment can be quite small. This presents a risk that the real signal could be of a similar size to other spurious signals. This App note discusses common sources of spurious variation in capacitance, and how to design your experiments to avoid them.

A large signal is always easier to interpret, so we recommend that the experiment is designed, and the sample selected, so that at least half of the cell's design range is used. Where that is not possible, extra care needs to be taken to ensure the data captured is of the best possible quality, as described in this document.

## Converting capacitance to force or displacement

The capacitance of a parallel plate capacitor is given by:

$$C = \frac{\epsilon_0 A}{d_{gap}}$$

Where  $C$  is the capacitance,  $\epsilon_0$  is the permittivity of free space (8.854 pF/m),  $A$  is the plate area and  $d$  the distance between them. We slightly modify this equation when drawing our calibration curves, so we use:

$$C = \frac{\alpha}{d + d_0} + C_p \qquad C = \frac{\alpha}{f + f_0} + C_p$$

For force and displacement sensors respectively. Here we have included a parallel capacitance term  $C_p$  to account for imperfections in

the capacitor, and we have replaced  $d_{gap}$  with  $d + d_0$ , where  $d_0$  is the gap at zero displacement and  $d$  is the change in gap that we want to measure. We have also combined  $\epsilon_0$  and  $A$  into a single coefficient  $\alpha$  for convenience. The second equation is similar, except for force sensors the first term is multiplied by a factor for the stiffness of the spring which converts the force applied to the sensor into a displacement for the capacitor to measure - but this too has been combined into  $\alpha$  for convenience. This does mean that for a displacement capacitor, the units of  $\alpha$  are given in  $\mu\text{m}\cdot\text{pF}$  whereas in the second equation it has units  $\text{N}\cdot\text{pF}$ .

The curves and equations provided on the first page of the calibration sheets for each sensor are obtained by fitting the equations above to measured points. These measurements are made at room temperature and in air, but the values can change due to various environmental factors. The rest of this document includes a description of the main factors, and on page 13 there is our recommended method for correcting for them.

All cells are supplied with a calibration curve for their primary response (i.e. displacement or force) and a temperature calibration<sup>1</sup>. All calibration data is archived at Razorbill Instruments, so if the factory calibration has been lost please contact razorbill instruments with the model type and serial number of your cell and we will be happy to supply a copy.

## Methods for on-site (re)calibration

The method of calibration depends on the type of cell. Force and distance calibrations are provided by Razorbill Instruments, and cells can be returned to the factory if updated calibrations are required. If you wish to perform your own calibration for force or displacement, please get in touch with Razorbill Instruments for more information about the parts and equipment required.

Calibrating a cell for temperature dependence requires no extra parts and is strongly recommended, even for cells that come with a temperature calibration from the factory. On-site calibration allows

---

<sup>1</sup> This applies to all cells sold in 2021 or later, before that only some models came with a temperature calibration.

variation in measurement equipment, cables, or thermometry to be considered.

Displacement sensing cells should be cooled and rewarmed with the supplied titanium calibration sample attached to prevent any movement of the displacement sensor. It is also helpful to “zero” the cell by applying 120V to both the tension and compression stacks, then reducing to 0V over about 30s. This gives the stacks the same history, and so the most similar thermal expansion behaviour, and reduces the force on the stiff sample.

Force-measuring cells can be temperature-calibrated simply by cooling and rewarming the cell with no sample attached, as that guarantees zero force. If your cell is fitted with both displacement and force sensors, two calibration experiments will be needed to calibrate each sensor separately.


Before taking a measurement, the cell should be held at a fixed temperature for some time to allow all the components in the cell to reach the desired temperature. It is recommended that measurements are taken both on cooling and warming, if hysteresis is observed the hold period was too short. Two hours is suggested as an initial hold period, but this may depend significantly on the design of your cryostat.

## FACTORS AFFECTING THE CALIBRATION

There are several factors affecting the calibration of the cells, some of which can be corrected for by adjusting the coefficients in the equations on page 3, and others which cannot be easily corrected so must be avoided instead.

### Temperature

Temperature has several affects on capacitance. It causes the plates to shrink, and also affects the gap between and alignment of the plates.



**Figure 1: An example temperature calibration curve. Total capacitance variation is c. 10 %, which corresponds to a significant spurious force reading. Note that the left axis does not go all the way to 0 as there is no correction for humidity. The left axis is the apparent force measured at each temperature in vacuum, converted from capacitance using the calibration curve measured at room temperature in air.**

When we consider a typical temperature calibration, we find that the change in capacitance between room temperature and 4 K is typically around 3-10 %. This is much greater than the 0.3 % attributable to change in plate area, so must derive mostly from plate movement. Razorbill Instruments has conducted several experiments to calibrate both force and displacement sensors at low temperature, and we have observed that the calibration curves remain fairly accurate if the change with temperature is treated as an offset in the quantity being measured. This means converting the capacitance measured during temperature calibration to a force or displacement, and subtracting that from the force or displacement measured.

For force sensors, there is an additional contribution which affects the gain of the sensor as well. The sensor works by using a stiff titanium beam to convert the force to a displacement, then measures that displacement. As the titanium gets stiffer on cooling, this changes the

gain of the sensor. Razorbill has measured this stiffening on a number of force sensing cells, and the change in  $\alpha$  is:

$$\alpha(T) = \alpha_{290K} \times (0.91 + 5 \cdot 10^{-5}T + 9 \cdot 10^{-7}T^2)$$


## Rate of change of temperature

Rapid changes in temperature lead to erroneous changes in capacitance. This occurs because different parts of the capacitor are at different temperatures. Let us imagine the at a cell is at 300 K and is cooled. Firstly, the titanium body of the cell cools, causing it to contract, and increasing capacitance. Next, insulating layers behind the capacitor plates cool, opening the gap, and decreasing capacitance. Finally, the capacitor plates cool and contract in area and thickness causing further decrease in capacitance. The actual process is more complex, involving complex thermal gradients as heat flows through the cell but the concept of uneven heating or cooling perturbing the capacitance remains valid. The extent of this effect depends on heat flow within the cryostat and through the cell which will vary significantly according to the presence or absence of exchange gas, the radiation environment, and the way heat is conducted into the cell.

The effects of rate of change of temperature on capacitance cannot easily be predicted and compensated for. It is therefore recommended that important measurements are only taken when the temperature is constant, and the cell has had time to reach thermal equilibrium.

Where it is necessary to take measurements during temperature changes (for example, when controlling the drive voltage to keep the sample at zero force during cooling) it is important to change temperature slowly enough that the measurement remains accurate. Practically we find it is best to change the temperature in steps of e.g. 30K, and then wait for thermal equilibrium and adjust the voltage at each step. Otherwise, if the (spurious) rate effect exceeds the genuine change due to differential expansion, the control loop will do more harm than good.






**Figure 2: A data sample from our QA process showing capacitance transients caused by rapid temperature changes (2 hour hold periods with 5 K/min ramps). This cell design (CS200T) has a lower temperature dependence than some others, making the transients more visible.**

Without intimate knowledge of each and every customer cryostat, we cannot easily recommend a maximum rate of change of temperature. We recommend cooling and rewarming the cell at your desired rate while monitoring hysteresis in the capacitance. If the hysteresis is unacceptably large, the rate must be reduced.

## Exposure to atmosphere

The capacitors show some dependence on atmospheric conditions, believed to be related to the adsorption of moisture from the air. When transferred into the cryostat, the capacitors become more stable. There is also a small change as the permittivity of air is slightly higher than that of vacuum.






**Figure 3: temperature-capacitance hysteresis for a cell in our QA cryostat. Hysteresis of c. 1 % during 5 K/min temperature ramps causes the blue (falling temperature) and red (rising temperature) curves to separate. Holding each temperature step for 2 hours causes rising and falling curves to meet. This data is used to generate the temperature calibration curves such as the one in Figure 1.**

When the cell is exposed to vacuum, the capacitance will drop slightly. This drop may take some time to stabilise fully. We recommend that 24 hours is allowed between pumping out the cryostat and starting



experiments. To minimise this effect, it is recommended that the cells are stored in a dry environment in their original cases when not in use. It may also be worth baking the cells at 60-90 °C for an hour before loading into the cryostat, if the sample allows this. Storing the cells in a vacuum chamber or dessicator can also help desorb moisture.



**Figure 4: repeated pump-downs on the same cell, starting from ambient each time. The capacitance initially varies according to atmospheric conditions, but after 24 - 48 hours in vacuum the variation is much reduced.**

If the cryostat vacuum is poor, or exchange gas is contaminated with air, water ice may form on the cell. Even very small quantities forming in the capacitor will greatly affect the measurement. If spikes in capacitance or other unusual behaviour are observed around 273 K, check the cryostat for leaks and confirm that good vacuum has been maintained. Similar issues can occur if other gasses, such as nitrogen, are allowed to condense in the cryostat, though the effect is smaller as their dielectric constants are much smaller than liquid water.

## Hysteresis

The displacement sensors show no detectable hysteresis. It should be remembered however that piezoelectric stacks do show hysteresis, so if the cell is held at constant drive voltage some creep of the measured value will be observed. This corresponds to genuine change in sample strain.

The Force sensors show some slight hysteresis, generally less than 0.5% of their full scale range<sup>2</sup>. This is worse when the cell is taken to its maximum load than when used at lower force. Consult the individual product datasheet for more information.

## Effect of displacement on a force sensor

Due to the miniaturised nature of the cells, and their flexure-guided construction, the operation of the piezoelectric stacks can cause a small distortion in the force-sensing area of the cell. If a force-sensing cell is operated at its maximum stroke, a small change in capacitance may be observed, but this is limited to about 0.1 %<sup>3</sup> of the full-scale range.

In a real experiment, the stiffness of the sample is such that the real signal is much larger, so this is rarely an issue except in cases where an extremely soft sample is used.

## Effect of force on a displacement sensor

The displacement sensor is positioned as close as reasonably practical to the sample but does not directly measure the extension of the sample itself. Extension of the sample mounting plates and epoxy are included in the measurement, and can be as large as the extension of the sample. The force transmitted through the cell's body will also cause some distortion, and in extreme cases of overloading the capacitor may be forced away from parallel. Without knowing the position and geometry of the sample, it is difficult to estimate the impact of this effect on the measurement. You can read more about the limitations of displacement sensing compared to force sensing in app note AP007.

The loss in accuracy due to these effects can be significant, so where possible it is advisable to calibrate the strain in the sample with an external measurement, such as the strain visible in a microscope image. Where external strain measurements are not possible, the following steps

---

<sup>2</sup> Applies to all cells except some early model FC100s, where it may be up to 3% FSR

<sup>3</sup> Older FC100 cells may be up to 1 %

will help minimise the difference between sample  $\Delta L$  and the displacement reported by the cell.

- The sample should be as long (and thus soft) as practicable
- The sample cross section should be small (and thus soft) relative to the sample plates or mounting arrangements
- The force should be kept within the maximum rating given in the cell datasheet, ideally below one third of the rating.

## Other sources of capacitance change

The FC100 shows some dependence on the tightness of the sample mounting screw. For more information refer to the FC100 datasheet. Other cells are not affected in this way.

## ADVICE ON EXPERIMENTS

To keep the effects mentioned above as small as possible, and to get the best achievable accuracy, we recommend you:

- Select your sample geometry to use at least half of the cell's design range, if possible
- Keep the sample stiffness within the design envelope of the cell. For displacement-sensing cells, lower is better.
- Ensure a good vacuum is formed in your cryostat, or, if using exchange gas, that the gas is not contaminated with air.
- Allow for a period of stabilisation after the cryostat is pumped out, 24 hours are recommended where practical.
- Allow for a period of stabilisation after a temperature change, 2 hours are recommended, depending on cryostat design.
- Keep the rate of temperature change under 4 K/min, or 0.5 K/min if the cell is driven in a closed loop configuration.

## Finding the location of zero

Even with the careful application of calibration data, it can be difficult to know the capacitance that corresponds to zero sample strain.

On a displacement sensing cell, differential thermal expansion between the Titanium calibration specimen, actual specimen, and epoxy will disturb the zero point. Best practice is to compare the measurements taken on a strained sample with those taken on an unstrained specimen.

On a force sensing cell, it is sometimes helpful to deliberately break the specimen at the end of the experiment (by applying the maximum tension voltage). This allows a clear no-load condition to be recorded with the minimum possible opportunity for error.

## Converting capacitance to displacement or force, correcting for temperature and other factors.

Taking all of the information in the above sections into account, the recommended procedure for converting the measured capacitance to force or displacement is:


1. Where an accurate zero point is required, consider baking out the cell (max 90°C) or storing it in a vacuum desiccator for 48 hours before loading into the cryostat.
2. Load the cell into the cryostat and pump for as long as practical, Razorbill suggests 24 hours.
3. Note the capacitance and convert it to force or displacement using the equation on the calibration curve. This will be your zero-point offset for this experiment.
4. When conducting your experiment, record the capacitance and convert it to force or displacement by:
  - a. Using the formula on the calibration sheet. If it is a force sensor, use the temperature dependent  $\alpha$  value described above.
  - b. Look up the change in  $d_0$  or  $f_0$  between the temperature you are at and the temperature at which you completed step 3. If your calibration curve does not have  $d_0$  or  $f_0$  (curves produced before 2023), convert from capacitance using the equation on the first page of the sheet.
  - c. Subtract the change from 4b and the value from 3 from the value calculated in 4a.

# TROUBLESHOOTING

The following symptoms occur from time to time, and can sometimes be easily solved:

## Small steps in capacitance (less than 5 %)

Small steps that are not associated with real sample behaviour are usually caused by changes in grounding arrangements. As it is common to use low thermal conductivity cable in cryostats, the cable braid and capacitor guards are often separated from ground by 3-10  $\Omega$ . If the coax connectors of the cell touch each other or the cryostat body, the AC potential of the braids and shields may change. In some configurations, parasitic capacitances between cable core and braid, braid to braid, and braid to core can appear in parallel with the measured capacitance and cause an offset. This offset is usually small enough to be ignored, as long as it is stable. But changes in grounding cause it to change.



**Figure 5: Steps in capacitance caused by intermittent changes in grounding. In this case the cable braid was going open circuit at the MMCX connector attached to the cell.**

Potential solutions include insulating over the coax connectors with Kapton or Teflon tape, or, if the electrical noise levels are acceptable, wiring the connectors to the cryostat metalwork.

The above-mentioned effect is exacerbated if one of the braids loses its connection to ground. Cables should be carefully inspected where

they join the connectors, and if frayed or appearing loose then the connector should be replaced. If a connector is suspected of poor mating, it is helpful to wrap it in uninsulated copper wire. If this resolves the issue, the connector and/or its partner should be replaced.

## Measured capacitance falls to zero

It is helpful to check the excitation voltage. If this has also collapsed close to zero a short circuit between the excitation (high) and ground or excitation and sense (low) has occurred. A multimeter may assist in finding the short circuit, but it's best to avoid directly probing the MMCX connectors as they are easily damaged. Check for stray conductive material in the connectors. A short circuit can be caused by foreign material entering the capacitor, and it is sometimes possible to clean it out, contact Razorbill Instruments for advice.

If the excitation voltage remains at its normal level, there is probably an open connection at some point in the cryostat wiring or inside the cell. Inspect the connectors for damaged centre pins – damage occurs easily if a probe or other object has been inserted, and check the cryostat wiring for continuity. If that doesn't resolve the problem, contact Razorbill Instruments for advice.

## Capacitor doesn't respond to changes in piezo voltage

This can have several causes:

- On a force-sensing cell, the sample is broken or detached
- Moving parts of the cell have been accidentally glued together
- The piezos are damaged or unplugged
- The same voltage in the same polarity is applied to all the piezoelectric stacks – the cells work on differential voltage.
- A glue joint within the cell is broken

Mounting an easy to handle sample (such as a piece of titanium foil) can help track down the issue.

## Hysteresis

Hysteresis with respect to temperature, if not originating from actual specimen behaviour, is best addressed by reducing the rate of change of temperature or increasing the hold time before a measurement is taken.

Hysteresis of displacement or force with respect to applied voltage is usually genuine and originates in the piezoelectric stacks. Treat the displacement or force as accurate and adjust the applied voltage accordingly.

Hysteresis of sample properties with respect to force or displacement should be small, as long as the sample is within the elastic limit of the material. If it is larger than expected, it may be a property of the sample or the epoxy used to secure it.

## Locating a fault

Before attempting a repair or returning a cell to Razorbill it can be helpful to confirm if the fault is in the wiring or in the cell itself. This can be done by substituting the cell for a dummy load. A 1 pF NPO or COG capacitor is a suitable load for cryogenic use. Additionally, the cell may be tested outside the cryostat and connected directly to your capacitance meter.

## Helping us to help you

If you are struggling to get good measurements, the team at Razorbill Instruments can give advice. It helps us to have as much information as possible. In particular, please try and capture the problematic behaviour on a datalogger or computer. Alongside capacitance it is helpful to record the quadrature (loss) component of the measurement and the excitation voltage.



It is also helpful to measure the resistance between components with a multimeter – ideally one that can measure impedances in the MΩ and higher. With the cables disconnected from the capacitance meter measure resistance for all positions in the following grid;

	Cryostat ground	High coax core	High coax shield	Low coax core
Low coax shield				
Low coax core				
High coax shield				
High coax core				

All these measurements should be open circuit. Please measure on a SMA or BNC connector, as pressing a multimeter probe against the middle pin of the small MMCX connector can permanently damage the connector.

