



Sample contacts.
Four 1x0.35mm



Cable. Copper twisted pair (pre Feb 2020) or PEEK woven loom (Feb 2020 onwards)
See datasheet for details and connector




Optional thermometer or low power heater below cable entry

Long sample

Typical Sample

Short sample




When very short samples are used, the WP100 must be mounted at an angle and the sample contacts will be slightly further from the sample.

This drawing applies to version 3 of the WP100. These have black PCBs. For drawings of earlier models (with blue or white PCBs) contact Razorbill Instruments.

Tolerances 0.5mm unless otherwise indicated

Drawing shows WP100L. WP100R is mirror image.



THIRD ANGLE

RAZORBILL INSTRUMENTS

PHONE: 0131 6688 488 EMAIL: jack@razorbillinstruments.com
UNLESS NOTED OTHERWISE, ALL DESIGNS AND DRAWINGS ARE CONFIDENTIAL

TITLE:

WP100 Customer Dwg

VERSION: 3 DATE: 2020-04-22

DWG NO. **WP100_cstmr_v3**

MATERIAL: Mixed

SHEET 1 OF 1

SCALE:5:1 PRINT ON A4
DO NOT SCALE DRAWING