


Sample contacts.
Four 1x0.35mm

Cable. Copper twisted
pair in PEEK woven loom
See datasheet for length
and connector

M2 x5 countersunk screw (included)

Tolerances 0.5mm unless
otherwise indicated

Drawing shows WP101L.
WP101R is mirror image.

 THIRD ANGLE		RAZORBILL INSTRUMENTS <small>PHONE: 0131 258 4329 EMAIL: info@razorbillinstruments.com</small> <small>UNLESS NOTED OTHERWISE, ALL DESIGNS AND DRAWINGS ARE CONFIDENTIAL</small>	
TITLE: <h1 style="margin: 0;">WP101 Customer Dwg</h1>			
VERSION: 1	DATE: 2022-12-01	DWG NO. WP101_cstmr	
MATERIAL: Mixed		SHEET 1 OF 1	SCALE:5:1 PRINT ON A4 DO NOT SCALE DRAWING